




**BELMONT
MEDICAL
TECHNOLOGIES**

Saving Lives. Together.

ALLON[®]

Safe, Precise, and Effective^{1,2}

PATIENT TEMPERATURE MANAGEMENT



CONTROL AND MAINTAIN NORMOTHERMIA THROUGHOUT THE ENTIRE PERIOPERATIVE CONTINUUM



PRE-OPERATIVE

Initiate patient warming before the procedure



INTRA-OPERATIVE

Manage patient temperature during the procedure



POST-OPERATIVE

Maintain patient body temperature in the recovery room

ALLON[®] ADVANCED TEMPERATURE MANAGEMENT

Allon is a non-invasive targeted temperature management system that controls and maintains patient body temperature throughout all stages of the perioperative continuum.

The system consists of an algorithm driven heat pump that supplies water through a closed-loop system to the ThermoWrap[®] disposable patient garment. Core and surface temperature sensors provide continuous patient temperature feedback to enable precise and automatic thermal regulation and monitoring.




THERMOWRAP[®]

ThermoWrap is a single use, one-piece flexible garment that acts in conjunction with the Allon control unit. ThermoWrap's flexible configurations facilitate the wrapping of individual parts of the body to maximize surface coverage without interfering with surgical sites. The CardiacWrap (shown next page, upper right) offers full access to the torso, making it a popular choice for cardiac procedures.


ThermoWrap is designed to conveniently attach to the control unit using two Quick Coupling Connectors— one for inflow and one for outflow. Integrated channels within the garment circulate water in a closed-loop system over a large surface area. ThermoWrap is secured via pressure sensitive adhesive strips which adhere to the garment, and not the patient.

YOUR PATIENT'S FAST-TRACK TO THE PACU^{3,4}

- Precise and automatic temperature control using proprietary algorithms
- Full-color, touch-screen display
- Simple, easy to use interface with intuitive controls
- Real-time monitoring and graphical display of temperature measurements
- Quick and easy set up
- Track and externally access procedure data using the CliniLogger™



Soft mesh garment with built-in pressure relief channels



DESIGNED WITH YOUR PATIENT IN MIND

- One-piece body shaped garment with various sizes and configurations
- Up to 85% body surface area coverage
- Soft mesh that is easy to wrap and secure
- Clutter-free with only one pair of connecting tubes

The ThermoWrap is available in a range of sizes and designs to accommodate your surgical needs and patient's anatomy:

Wrap Type	Part Number	Quantity Per Box	Patient Height/Weight	Wrap Length/Width (m)
Cardiac	512-03363	12	165-182 cm Fits all adults	1.348 / 1.319
Adult/Universal	512-03166	12	168-180 cm (over 66")	1.904 / 1.321
	512-03160	12	152-168 cm (60-66")	1.934 / 1.295
	512-03153	12	135-152 cm (53-60")	1.744 / 1.212
Pediatric	512-03148	12	122-135 cm (48-53")	1.582 / 1.193
	512-03141	12	104-122 cm (41-48")	1.398 / 1.068
	512-03136	12	91-104 cm (36-41")	1.225 / 0.841
Infant	524-03125	24	7-11 kg (15.5-24.3 lb)	0.983 / 0.629
	524-03121	24	4-7 kg (8.8-15.5 lb)	0.698 / 0.604
	524-03118	24	2.5 to 4 kg (5.5-8.8 lb)	0.660 / 0.465

1 Kjellman et al.: Comparing ambient, air-convection, and fluid-convection heating techniques in treating hypothermic burn patients, a clinical RCT. Annals of Surgical Innovation and Research 2011 5:4.
 2 Motta et al.: Effectiveness of a Circulating-Water Warming Garment in Rewarming After Pediatric Cardiac Surgery Using Hypothermic Cardiopulmonary Bypass. J Cardiothoracic and Vascular Anesthesia 2004. 18:2, 148-151.
 3 Rohrer B., Penick E., Zahedi F., and et al., Comparison of forced-air and water circulating warming for prevention of hypothermia during transcatheter aortic valve replacement. PLoS one, 2017; 12(6): e0178600.
 4 Neshar, N., Uretzky, G., Insler, S., Nataf, P., Frolkis, I., Pineau, E., ... Weinbroum, A. A. (2005). Thermo-wrap technology preserves normothermia better than routine thermal care in patients undergoing off-pump coronary artery bypass and is associated with lower immune response and lesser myocardial damage. The Journal of Thoracic and Cardiovascular Surgery, 129(6), 1371-1378. doi:10.1016/j.jtcvs.2004.08.021

ORDERING INFORMATION

CONTROL UNIT

Allon® 2001 (115V)	Part #100-00001
Allon® 2001 (230V)	Part #100-00002

ACCESSORY KIT

(Select one only for each Allon unit)

ADULT/INFANT KIT FOR DISPOSABLE SENSORS	Part #200-00410
--	------------------------

The Kit contains:

• Core Adapter Cable	Part #014-00028
• Surface Adapter Cable	Part #014-00129
• Connecting Water Tubes (2 by 2 way)	Part #200-00109

REUSABLE KITS:

(NOT AVAILABLE IN THE UNITED STATES AND SELECT MARKETS)

ADULT KIT FOR REUSABLE SENSORS*	Part #200-00400
--	------------------------

The Kit contains:

• Reusable Adult Core Temperature Probe	Part #014-00020
• Reusable Surface Temperature Probe	Part #014-00021
• Filter unit	Part #200-R0130
• Connecting Water Tubes (2 by 2 way)	Part #200-00109

INFANT KIT FOR REUSABLE SENSORS*	Part #200-00420
---	------------------------

The Kit contains:

• Reusable Infant Core Temperature Probe	Part #014-00005
• Reusable Surface Temperature Probe	Part #014-00021
• Filter unit	Part #200-R0130
• Connecting Water Tubes (2 by 2 way)	Part #200-00109

ACCESSORIES

Disposable Core Temperature Probe	Part #014-00322
Disposable Surface Temperature Probe	Part #014-00321
Core Adapter Cable	Part #014-00028
Surface Adapter Cable	Part #014-00129
Reusable Core Temperature Probe, Adult*	Part #014-00020
Reusable Core Temperature Probe, Infant*	Part #014-00005
Reusable Surface Temperature Probe*	Part #014-00021
CliniLogger™	Part #017-00250

*Reusable temperature probes are not for sale in the United States.

HARDWARE SPECIFICATIONS

HEAT EXCHANGERS

Peltier Technology (TECs)

WATER TANK

Sterile or 0.22 micron filtered tap water
Tank Capacity: 6 liters (1.6 gal.)

WATER TEMPERATURE

Water Temperature Accuracy: $\pm 0.3^{\circ}\text{C}$ (0.4°F)
Water Temperature (outflow): $13\text{--}40.8^{\circ}\text{C}$ ($55.4\text{--}105.4^{\circ}\text{F}$)

PUMP

Water Circulating Pump
Pump Rate: 0.2–1.25 L/min
Protected by 263 micron filter

PATIENT TEMPERATURE CHANNELS

2 Channels: Core, Surface
Sensor Temperature Range: 15°C to 44°C (59°F to 111.2°F)
Sensor Temperature Accuracy: $\pm 0.3^{\circ}\text{C}$ (0.4°F)

TEMPERATURE AND PRESSURE SENSORS

System Sensors:

- 3 Internal Temperature Sensors: Water in, Water out, Thermostat
- 2 Pressure sensors

Safety Measures:

- Over pressure protection and alarm
- High water temperature protection and alarm

PHYSICAL DIMENSIONS

Mobile unit with 4 wheels and 2 brakes
260 mm W x 625 mm D x 940 mm H /
(10.23"W x 24.6"D x 37"H)

NET WEIGHT

34 kg / 75 lb

INPUT POWER

115VAC 5.7 A 50/60 HZ
230VAC 2.9 A 50/60 HZ

THERMOWRAP

Duration of use: up to 28 hours

LCD DISPLAY

Size: 5.7" color display
Resolution: 320x240

For more information about the products we offer,
contact us at sales@belmontmedtech.com.



PATIENT
WARMING



PATIENT
COOLING



FLUID
MANAGEMENT



780 BOSTON ROAD
BILLERICA, MA 01821, USA
USA: 866.663.0212
WORLDWIDE: +1 978.663.0212
FAX: 978.663.0214
WWW.BELMONTMEDTECH.COM

