



**ALLON<sup>®</sup>**

Safe, Precise, and Effective<sup>1,2</sup>

PATIENT TEMPERATURE MANAGEMENT





CONTROL  
AND MAINTAIN  
NORMOTHERMIA  
THROUGHOUT  
THE ENTIRE  
PERIOPERATIVE  
CONTINUUM



## PRE-OPERATIVE

Initiate patient warming  
before the procedure



## INTRA-OPERATIVE

Manage patient temperature during  
the procedure



## POST-OPERATIVE

Maintain patient body temperature in  
the recovery room

# ALLON® ADVANCED TEMPERATURE MANAGEMENT

Allon is a non-invasive targeted temperature management system that controls and maintains patient body temperature throughout all stages of the perioperative continuum – without impacting laminar air flow and potentially increasing the risk for surgical site infections.

The system consists of an algorithm driven heat pump that supplies water through a closed-loop system to the ThermoWrap® disposable patient garment.



## THERMOWRAP®

ThermoWrap is a single use, one-piece flexible garment that is used in conjunction with the Allon control unit. ThermoWrap's flexible configurations facilitate the wrapping of individual parts of the body to maximize surface coverage without interfering with surgical sites. The CardiacWrap (shown next page, upper right) offers full access to the torso, making it a popular choice for cardiac procedures.

## PREDICTABLE WARMING – ALL OF THE TIME



Not every patient or procedure is at the same risk for a surgical site infection (SSI). Longer surgical procedures, transplant, orthopedic and certain cardiac surgeries all have significantly higher risks for SSIs compared to other shorter, less complex surgeries.<sup>5</sup>

When forced air warming devices are used, the laminar airflow is impacted, and doubt remains whether the ultra clean air ventilation systems are able to properly restore air flow. As a result, bacteria can be blown into the surgical site, potentially causing an infection.

With a reputation for predictable warming, Allon has been proven to maintain normothermia even in complex cases. A conductive warming device, Allon does not disturb laminar air flow, nor increase the odds of a surgical site infection which is crucial – especially for high risk patients and populations.

## YOUR PATIENT'S FAST-TRACK TO THE PACU<sup>3, 4</sup>

- Precise and automatic temperature control using proprietary algorithms
- Full-color, touch-screen display
- Simple, easy to use interface with intuitive controls
- Real-time monitoring and graphical display of temperature measurements
- Quick and easy set up
- Track and externally access procedure data using the CliniLogger™



## DESIGNED WITH YOUR PATIENT IN MIND

- One-piece body shaped garment with various sizes and configurations
- ThermoWrap provides up to 80% surface area coverage\*
- Soft mesh that is easy to wrap and secure
- Clutter-free with only one pair of connecting tubes

\* Universal ThermoWrap only.

ThermoWrap is available in a range of sizes and designs to accommodate your surgical needs and patient's anatomy:

Wrap Type	Part Number	Quantity Per Box	Patient Height/Weight	Wrap Length/Width (m)
Cardiac	512-03363	12	165-182 cm Fits all adults	1.348 / 1.319
Adult/ Universal	512-03166	12	168-180 cm (over 66")	1.904 / 1.321
	512-03160	12	152-168 cm (60-66")	1.934 / 1.295
	512-03153	12	135-152 cm (53-60")	1.744 / 1.212
Pediatric	512-03148	12	122-135 cm (48-53")	1.582 / 1.193
	512-03141	12	104-122 cm (41-48")	1.398 / 1.068
	512-03136	12	91-104 cm (36-41")	1.225 / 0.841
Infant	512-03131	12	79-91 cm (31-36")	1.118 / 0.739
	524-03125	24	7-11 kg (15.5-24.3 lb)	0.983 / 0.629
	524-03121	24	4-7 kg (8.8-15.5 lb)	0.698 / 0.604
	524-03118	24	2.5 to 4 kg (5.5-8.8 lb)	0.660 / 0.465

1 Kjellman et al.: Comparing ambient, air-convection, and fluid-convection heating techniques in treating hypothermic burn patients, a clinical RCT. Annals of Surgical Innovation and Research 2011 5:4.  
2 Motta et al.: Effectiveness of a Circulating-Water Warming Garment in Rewarming After Pediatric Cardiac Surgery Using Hypothermic Cardiopulmonary Bypass. J Cardiothoracic and Vascular Anesthesia 2004. 18:2, 148-151.  
3 Rohrer B., Penick E., Zahedi F., and et al.: Comparison of forced-air and water circulating warming for prevention of hypothermia during transcatheter aortic valve replacement. PLoS one, 2017; 12(6): e0178600.  
4 Neshor, N., Uretsky, G., Insler, S., Nataf, P., Frolkis, I., Pineau, E., ... Weinbroum, A. A. (2005). Thermo-wrap technology preserves normothermia better than routine thermal care in patients undergoing off-pump coronary artery bypass and is associated with lower immune response and lesser myocardial damage. The Journal of Thoracic and Cardiovascular Surgery, 129(6), 1371-1378. doi:10.1016/j.jtcvs.2004.08.021  
5 Mu, Yi et al. Improving Risk-Adjusted Measures of Surgical Site Infection for the National Healthcare Safety Network. Infection Control and Hospital Epidemiology, Vol. 32, No. 10 (October 2011), pp. 970-986.



## ORDERING INFORMATION CONTROL UNIT

Allon® 2001 (115V)  
Allon® 2001 (230V)

Part #100-00001  
Part #100-00002

## ACCESSORY KIT

(Select one only for each Allon unit)

### ADULT/INFANT KIT FOR DISPOSABLE SENSORS

Part #200-00410

The Kit contains:

- Core Adapter Cable Part #014-00028
- Surface Adapter Cable Part #014-00129
- Connecting Water Tubes (2 by 2 way) Part #200-00109

### REUSABLE KITS:

(NOT AVAILABLE IN THE UNITED STATES AND SELECT MARKETS)

### ADULT KIT FOR REUSABLE SENSORS\*

Part #200-00400

The Kit contains:

- Reusable Adult Core Temperature Probe Part #014-00020
- Reusable Surface Temperature Probe Part #014-00021
- Filter unit Part #200-R0130
- Connecting Water Tubes (2 by 2 way) Part #200-00109

### INFANT KIT FOR REUSABLE SENSORS\*

Part #200-00420

The Kit contains:

- Reusable Infant Core Temperature Probe Part #014-00005
- Reusable Surface Temperature Probe Part #014-00021
- Filter unit Part #200-R0130
- Connecting Water Tubes (2 by 2 way) Part #200-00109

## ACCESSORIES

Disposable Core Temperature Probe, 9FR, Smiths Medical	Part #014-00035
Disposable Core Temperature Probe, 7FR, Metko	Part #014-00036
Disposable Core Temperature Probe, 9 FR, TE	Part #014-00038
Disposable Core Temperature Probe, 9FR, DeRoyal	Part #014-00220
Disposable Surface Temperature Probe	Part #014-00321
Core Adapter Cable	Part #014-00028
Surface Adapter Cable	Part #014-00129
Reusable Core Temperature Probe, Adult*	Part #014-00020
Reusable Core Temperature Probe, Infant*	Part #014-00005
Reusable Surface Temperature Probe*	Part #014-00021
CliniLogger™	Part #017-00250
Temperature Splitter Kit	Part #200-01200
The Temperature Splitter Kit contains:	
• One (1) Temperature Splitter	Part #017-01200
• Two (2) Temp-Out Adapter Cables	Part #014-00012

\*Reusable temperature probes are not for sale in the United States.

## HARDWARE SPECIFICATIONS

### HEAT EXCHANGERS

Peltier Technology (TECs)

### WATER TANK

Sterile or 0.22 micron filtered tap water  
Tank Capacity: 6 liters (1.6 gal.)

### WATER TEMPERATURE

Water Temperature Accuracy:  $\pm 0.3^{\circ}\text{C}$  ( $0.4^{\circ}\text{F}$ )  
Water Temperature (outflow):  $13\text{--}40.8^{\circ}\text{C}$  ( $55.4\text{--}105.4^{\circ}\text{F}$ )

### PUMP

Water Circulating Pump  
Pump Rate: 0.2–1.25 L/min  
Protected by 263 micron filter

### PATIENT TEMPERATURE CHANNELS

2 Channels: Core, Surface  
Sensor Temperature Range:  $15^{\circ}\text{C}$  to  $44^{\circ}\text{C}$  ( $59^{\circ}\text{F}$  to  $111.2^{\circ}\text{F}$ )  
Sensor Temperature Accuracy:  $\pm 0.3^{\circ}\text{C}$  ( $0.4^{\circ}\text{F}$ )

### TEMPERATURE AND PRESSURE SENSORS

System Sensors:

- 3 Internal Temperature Sensors: Water in, Water out, Thermostat
  - 2 Pressure sensors
- Safety Measures:
- Over pressure protection and alarm
  - High water temperature protection and alarm

### PHYSICAL DIMENSIONS

Mobile unit with 4 wheels and 2 brakes  
260 mm W x 625 mm D x 940 mm H /  
(10.23"W x 24.6"D x 37"H)

### NET WEIGHT

34 kg / 75 lb

### INPUT POWER

115VAC 5.7 A 50/60 HZ  
230VAC 2.9 A 50/60 HZ

### THERMOWRAP

Duration of use: up to 28 hours

### LCD DISPLAY

Size: 5.7" color display  
Resolution: 320x240

For more information about the products we offer,  
contact us at [sales@belmontmedtech.com](mailto:sales@belmontmedtech.com).



PATIENT  
WARMING



PATIENT  
COOLING



FLUID  
MANAGEMENT



780 BOSTON ROAD  
BILLERICA, MA 01821, USA  
USA: 866.663.0212  
WORLDWIDE: +1 978.663.0212  
[WWW.BELMONTMEDTECH.COM](http://WWW.BELMONTMEDTECH.COM)



Belmont Instrument, LLC  
Certified to ISO 13485:2016  
File Number: A18173

Belmont Medical Technologies™ is a trademark of Belmont Instrument, LLC  
2024 © Belmont Instrument, LLC | 701-00272 Rev. 001