

PEDIATRIC CUREWRAP®




For all Pediatric CureWraps® except Model 3500.
Pediatric CureWraps® require **2x2 Way** Connecting Water Hoses PN #200-00109.
(See Image A). For Model 3500, see Model 3500: Adult and Pediatric CureWrap®.

MODEL 3500: ADULT/PEDIATRIC CUREWRAP®

For Adult
CureWrap®
Model 3500
and Pediatric
CureWrap®
Model 3500.

For other
wrap models,
see Pediatric
CureWrap®.
Model 3500
CureWraps®
require **2x3
Way**
Connecting
Water Hoses
PN #200-00147
(See below).

Apply the CureWrap to the patient according to the instructions below and Pediatric CureWrap® or Model 3500 CureWrap® illustrations as appropriate.

- Unroll the CureWrap on the bed with the tubes towards the head of the bed with the lot and part numbers legible.
- Note the model number of the wrap being used to determine the appropriate wrap tubes. Pediatric CureWraps (all wrap models except 3500) require 2x2 Way Connecting Water Hoses PN #200-00109 (See Image A). Model 3500 for patients over 135 cm requires 2x3 Way Connecting Water Hoses PN#200-00147 as shown in Image #5. Model 3500 is sold in PN# 508-03500 in addition to PN# PED-XL008.
- Connect CureWrap to CritiCool by inserting the wrap tubes into the plastic ends of the blue Connecting Water Hoses, listening for the click. Confirm that the clamps on the wrap tubes are open. (See Image A and/or Image #5).
- Activate CritiCool.
- Position patient on CureWrap with the patient's head towards the tubes. Refer to Image B (and/or Images #1-5) for correct armpit positioning.
- Wait until CureWrap fills with water prior to securing it around the patient. Water will not flow without a valid core temperature. See CritiCool User Manual for instructions on temperature probes.
- When closing the CureWrap around the patient, ensure that it fits loosely.
- There should be a finger's width space between the patient and the CureWrap. The head and arms of infant patients should not be wrapped. For infant and pediatric patients, see Images C and D. For adult patients, refer to Images 6-10.
- Prevent any thermal isolation, such as a pillow or other items, between the Wrap and the patient's body.

CAUTION:

- U.S. Federal Law restricts this device for sale to or on the order of a physician.
- Avoid handling CureWrap with sharp objects.
- Do not use the CureWrap to lift patient.
- Ensure that CureWrap is not folded or clamped.
- If treatment exceeds 120 hours, replace CureWrap.

NOTE: In case of leakage, the CritiCool device should be turned off.

NOTE: Any serious incident that occurs in relation to this device should be reported to Belmont Medical Technologies and to the competent authority of the member state in which the user and/or patient is established.

To disconnect the CureWrap:

- Switch off CritiCool device or switch it to Standby Mode.
- Wait a few seconds.
- Close clamps and disconnect.
- Unwrap patient and dispose of CureWrap.
- See CritiCool user manual for a full list of warnings / contraindications.

Intended use - CureWrap is a single use device designed exclusively for use with the CritiCool system. It is designed to facilitate the wrapping of individual parts of the body (chest, arms, thighs, etc) to maximize surface coverage. CureWrap is intended to maintain preset body temperature as determined by a physician when used with the CritiCool system.

Patient Target Group – CureWrap can be used with infant, pediatric and adult patients.

Intended Users – Trained health care professionals, including nurses and physicians.

Contraindications - CureWrap should not come into contact with open wounds. Caution should be taken when using CureWrap with patients that have underlying skin conditions.

Residual Risk – Possible patient injury could occur if the garment is wrapped around the patient too tightly.

Device Disposal – CureWraps are designated as general medical waste. CureWrap should be disposed of in the general waste stream.

Risk of Reuse - CureWrap is for single use only. Reuse could lead to cross contamination and/or irritation.

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